

# PMUX-S (serial)

## Contents of package

- 1 PMUX-S 'high current cell multiplexer'
- 1 Controller card RMUX/P (plug-in card for IM6/6eX)
- 1 cable SUBD44 to SUBD25 (RMUX/P->PMUX)
- 1 IM connector cable set (banana->Lemosa, red/black, blue/green)
- 1 PP connector cable set (banana/banana red/black, banana/Lemosa blue/green)
- 16 sets of cell cables (4 lines each: red/black, blue/green)
- 1 RMUX cable set (Lemosa/Lemosa, only use if RMUX is used without PMUX)

## CAUTION

The PMUX-S is configured for sequentially measuring through up to 16 individual cells using either the internal potentiostat of the IM6/6eX **OR** an external potentiostat of the PP-series or an electronic load of the EL-series.



**Please take care that only one potentiostat EITHER the internal OR an external is connected to the PMUX-S. If both potentiostats are connected at the same time the system will get damaged.**




**Please take care that the maximum current is 5 A. Higher currents will damage the PMUX-S.**



**If you intend to use the RMUX without the PMUX-S the RMUX hardware setting must be changed. Please contact ZAHNER for more information. The RMUX cable set is used only for that application.**

## Hardware installation

Switch off the IM6 and remove the dummy front panels named 'IM6 extension' **RIGHT** to the 'AHSI-2' card. Push in the **RMUX** card in one of the empty extension slots and fix the screws. If the **EPC40** for connecting the PP is not installed, do it now in the same way.



### Connection of modules when using the internal IM-POT

1. Connect the PMUX-S to the RMUX/P using the multi-pin Sub-D/Sub-D cable. Connect the 44-pin (3-row) Sub-D connector to the RMUX/P and the 25-pin Sub-D connector to the CTRL terminal of the PMUX-S.
2. Connect the PMUX-S to the IM6/6eX using the two banana/Lemosa cables red/black and blue/green. Connect the banana plugs (red/black and blue/green) to the outlets  $RE_{IM}$  (green),  $TES_{IM}$  (blue),  $CE_{IM}$  (red), and  $TE_{IM}$  (black) of the PMUX-S. Connect the Lemosa plugs to the **Probe-I** and **Probe-E** outlets of the IM6/6eX.
3. **LEAVE** all connectors in the section **PP System** ( $TES_{IM}$ ,  $RE_{PP}$ ,  $CE_{PP}$ , and  $TE_{IM}$ ) **UNCONNECTED !!**
4. Connect the channels you need (1 – 16) using all four lines each (RE, TES, CE, and TE) to the according electrodes or points of your cells. Channel 1 to cell 1, channel 2 to cell 2, etc.

RE = reference electrode (green)  
TES = test electrode sense (blue)  
CE = counter electrode (red)  
TE = test electrode (working electrode) (black)

### Connection of modules when using an external POT

1. Connect the PMUX-S to the RMUX/P using the multi-pin Sub-D/Sub-D cable. Connect the 44-pin (3-row) Sub-D connector to the RMUX/P and the 25-pin Sub-D connector to the CTRL terminal of the PMUX.
2. Connect the PP device to the EPC40 card as described in the PP manual.
3. Connect the PMUX-S to the PP using the thick banana/banana cable red/black and the banana/Lemosa cable blue/green. Connect the blue and green banana plugs of the banana/Lemosa cable to the  $TES_{PP}$  and  $RE_{PP}$  outlets of the PMUX-S and the Lemosa plug to the **Sense** outlet of the PP device. Connect the small banana plugs of the red/black to the outlets  $CE_{PP}$  (red),  $TE_{PP}$  (black) of the PMUX-S and the other end to the **Power** outlets of the PP.
4. **LEAVE** all connectors  $RE_{IM}$  (green),  $TES_{IM}$  (blue),  $CE_{IM}$  (red), and  $TE_{IM}$  (black) in the section **IM System UNCONNECTED !!**
5. Connect the channels you need (1 – 16) with all four lines each (RE, TES, CE, and TE) to the according electrodes or points of your cells. Channel 1 to cell 1, channel 2 to cell 2, etc.
6. Be careful not to exceed the maximum current of +5 A.