

# Dark Field Photoelectron Emission Microscopy of Micron Scale Few Layer Graphene



## Introduction



We present dark field Photoelectron Emission Microscopy and energy filtered electron emission spectromicroscopy with high spatial and wave-vector resolution on few-layer epitaxial graphene on SiC(000-1) grown by furnace annealing.

Conventional electron spectroscopy methods are limited in providing simultaneous real and reciprocal or k-space information from small areas under laboratory conditions. Therefore, the characterization of materials with only micron scale sample homogeneity such as epitaxially grown graphene requires new instrumentation. Recent improvements in aberration compensated energy-filtered photoelectron emission microscopy (PEEM) can overcome the known limitations in both synchrotron and laboratory environments.




Energy-filtered, threshold PEEM is used to locate distinct zones of FLG graphene. In each region, selected by a field aperture, the k-space information is imaged using appropriate transfer optics. By selecting the photoelectron intensity at a given wave vector and using the inverse transfer optics, dark field PEEM gives a spatial distribution of the angular photoelectron emission. Only the combination of high lateral, high energy, high k-resolution and controlled switching between real space and k-space allows detailed understanding of micron size sample sites with 1-3 layers graphene.

## Instrumental Details




Advantage of the NanoESCA lens design is that the analysis area for k-space analysis can precisely defined when the instrument is operated in real space and vice versa. Both apertures CA and FA are operated independently from each other. In the back focal plane of the first lens L1 the angular image is formed. In this plane a contrast aperture CA is situated. This CA allows to define the angular acceptance of the instrument. A large CA is chosen for a maximum k acceptance a small CA for dark field imaging. After the transfer lens L2 the first intermediate image is formed. In the plane of the first intermediate image a field aperture (iris- or a novel micro-aperture) can be activated to precisely define the local analysis area on the sample for space k-space imaging. Since the k-space and real space acceptance are defined in the front part of the instrument, the actual switching process has no influence on the area and angular definition. Switching is done by using the projective lenses to shift the focal plane of the image plane into the entrance of the energy analyser. As a result the NanoESCA lens design allows controlled switching between real space and k-space.


## Secondary Electron Spectra and Work-Function Map



## High resolution $\mu$ -ARUPS (from 5.7 $\mu\text{m}$ area on single graphene site)



## Darkfield PEEM (k-space images at E-Ef = 21.0 eV, recorded from a 73 $\mu\text{m}$ FoV using He I Excitation)



## He II: Early results using He II excitation

