

CryoAdvance Configuration/Information/Dimensions

www.montanainstruments.com

+1.406.551.2796

MONTANA INSTRUMENTS° COLD SCIENCE MADE SIMPLE



info@montanainstruments.com



System/Product Description

Basic description of CryoAdvance:

- ϕ 50mm or ϕ 100mm ID radiation shield/sample space
- Bare platform with options for radiation shield/vacuum housings and sample mounting
- Mounts to table parallel or 45° (Inch or Metric)
- Air cooled compressor
- Approx weight: 29Kg


Α

Β


С

D


• See Options for more detail



1



1



1

D

REVISION SHEET 3 OF 4

D

